

INSTRUCTIONS for Phos-tag<sup>®</sup> BTL-105 (Purification of Phosphorylated Peptide using Phos-tag<sup>®</sup> Biotin & Streptavidin-Agarose)  
by the NARD Institute Ltd. & MANAC Incorporated Group (<http://www.Phos-tag.com>)

## Introduction

Phos-tag<sup>®</sup> (a dinuclear zinc(II) complex) acts as a novel phosphate capture molecule in aqueous solution at physiological pH (6 ~ 8). The anion selectivity indexes of Phos-tag<sup>®</sup> are as follows: 16,000 for phenyl phosphate dianion, 3 for SO<sub>4</sub><sup>2-</sup>, 1 for CH<sub>3</sub>COO<sup>-</sup>, 0.02 for Cl<sup>-</sup>, <0.02 for diphenyl phosphate monoanion and alkyl sulfate (ROSO<sub>3</sub><sup>-</sup>) monoanion. Phos-tag<sup>®</sup> BTL-105 is a biotin derivative attached with Phos-tag<sup>®</sup> molecules, which allows the fast and effective purification of phosphate dianions (*e.g.*, phosphorylated peptides, phosphorylated proteins, phosphorylated lipids, *etc.*) using streptavidin-bound agarose (*Sigma-Aldrich Fine Chemicals*). To ensure best performance and trouble-free operation, please read the following instructions before using Phos-tag<sup>®</sup> BTL-105.

## Product Specifications for Phos-tag<sup>®</sup> BTL-105

Mol. Wt.: 880.1  
Binding site: Phos-tag<sup>®</sup>/biotin = 1:1  
Form: Zn<sup>2+</sup>-unbound ligand  
Storage place: at +2 ~ +8 °C in the dark

## Additional Materials Required

The followings are the materials required for use with Phos-tag<sup>®</sup> BTL-105.

- 1) *Zinc(II) Acetate Solution*  
0.50 mM Zn(CH<sub>3</sub>COO)<sub>2</sub> containing 5.0 mM Tris-acetate (pH 7.4)
- 2) *Tris acetate Buffer*  
5.0 mM Tris-acetate solution (pH 7.4)
- 3) *Balancing Buffer (pH 7.4)*  
5.0 mM Tris-acetate containing 10 μM Zn(CH<sub>3</sub>COO)<sub>2</sub>
- 4) *Washing Buffer (pH 7.4)*  
5.0 mM Tris-acetate containing 0.50 M NaNO<sub>3</sub>
- 5) *Elution Buffer (optional)*  
1.0 mM NaH<sub>2</sub>PO<sub>4</sub>-NaOH (pH 7.0) containing 0.50 M NaNO<sub>3</sub>
- 6) Micropipette for 0.30 mL
- 7) Centrifugal filtration units  
Membrane pore size: 0.22 ~ 0.45 μm  
Sample volume: 0.40 ~ 0.50 mL
- 8) Centrifuge for 2,000×g
- 9) Streptavidin-agarose (*Sigma-Aldrich Fine Chemicals*)

## Sample Preparation




- 1) Competing anions (*e.g.*, inorganic phosphate, thiolate) and metal chelating agents (*e.g.*, EDTA) should be removed from the sample as much as possible.  
Note: The protease inhibitors except of thiol compounds are not effective on the purification. The large amount of surfactant (*e.g.*, SDS) lowers the efficiency of the purification.
- 2) The sample is dissolved in the 5.0 mM Tris-acetate buffer (pH 7.4, 0.30 mL) at room temperature.
- 3) The final concentration of phosphorylated compounds in the sample solution (0.3 mL) should be below 15 nmol/mL.

## Purification Procedure


- 1) Apply *ca.* 0.3 mL of Streptavidine Agarose (in suspension) in the sample reservoir of the centrifugal filtration unit.
- 2) Centrifuge the unit at 2,000×*g* for 15 sec and the filtrate is discarded.
- 3) Apply 0.30 mL of the 5.0 mM Tris-acetate buffer in the sample reservoir.  
Centrifuge the unit at 2,000×*g* for 15 sec and the filtrate is discarded.  
(The operations are repeated 5 time)
- 4) Apply 0.30 mL of 0.12 mM Phos-tag<sup>®</sup> BTL-105 containing 0.50 mM Zn(CH<sub>3</sub>COO)<sub>2</sub> and 5.0 mM Tris-acetate (pH 7.4) in the sample reservoir.
- 5) Equilibrate for 5 min.
- 6) Centrifuge the unit at 2,000×*g* for 15 sec and the filtrate is discarded.
- 7) Apply 0.30 mL of the 5.0 mM Tris-acetate buffer in the sample reservoir.  
Centrifuge the unit at 2,000×*g* for 15 sec and the filtrate is discarded.  
(The operations are repeated 5 time)
- 8) Apply 0.30 mL of 10 μM Zn(CH<sub>3</sub>COO)<sub>2</sub> containing 5.0 mM Tris-acetate (pH 7.4) in the sample reservoir.
- 9) Centrifuge the unit at 2,000×*g* for 15 sec and the filtrate is discarded.
- 10) Add 0.3 mL of the sample solution into the sample reservoir.
- 11) Equilibrate for 5 min in order to sufficiently bind phosphorylated compounds to Phos-tag<sup>®</sup>.
- 12) Centrifuge the unit at 2,000×*g* for 15 sec. The filtrate contains non-phosphorylated compounds.
- 13) Add 0.30 mL of the washing buffer into the sample reservoir. Centrifuge the unit at 2,000×*g* for 15 sec. The filtrate contains non-phosphorylated compounds.  
(The operations are repeated 3 time)
- 14) Add 0.30 mL of the elution buffer into the sample reservoir. Centrifuge the unit at 2,000×*g* for 15 sec and the filtrate contains phosphorylated compounds.  
(The operations are repeated 2~4 time)
- 15) Analyze the obtained solutions containing phosphorylated compounds.

# Purification Protocol for Phosphorylated Peptides


using Phos-tag<sup>®</sup> Biotin & Streptavidin-Agarose




**Phos-tag® BTL-105**



Sample solution in 5.0 mM Tris-CH<sub>3</sub>COOH (pH 7.4, 0.30 mL)



p60c-src peptide 521-533 (p60c)  
14 nmol/mL



Phosphorylated p60c-src peptide 521-533 (P-p60c)  
12 nmol/mL

1) Separation result in the presence of Phos-tag® Biotin BTL-105

| Fraction No. | 1  | 2  | 3 | 4 | 5  | 6  |
|--------------|----|----|---|---|----|----|
| p60c (%)     | 67 | 24 | 9 | 0 | 0  | 0  |
| P-p60c (%)   | 0  | 0  | 0 | 0 | 83 | 17 |

← Washing → Eluting

2) Separation result in the absence of Phos-tag® Biotin BTL-105

| Fraction No. | 1  | 2  | 3 | 4 | 5 | 6 |
|--------------|----|----|---|---|---|---|
| p60c (%)     | 75 | 20 | 5 | 0 | 0 | 0 |
| P-p60c (%)   | 83 | 17 | 0 | 0 | 0 | 0 |

← Washing → Eluting

Quantitative analysis of the peptides is conducted by HPLC as follows.

Column: Shiseido Capcell Pak C<sub>18</sub> type UG80, 150 mm in length, 4.6 mm in  $\phi$ .

Eluent: 0.1%(v/v) trifluoroacetic acid, 14%(v/v) acetonitrile aqueous solution.

Detector: UV at 266 nm. Flow speed: 1 mL/min. Column temperature: 40 °C.

Retention time: 4.6 min for P-p60c; 12.6 min for p60c

## Surface Plasmon Analysis of Phosphorylated $\beta$ -Casein Using a Phos-tag<sup>®</sup> Biotin- and Streptavidin-Bound Sensor Chip

### Ligand binding procedures

- 1) A sensor chip, onto which streptavidin was bound (Sensor Chip SA, made by Biacore), was set in Biacore J (made by Biacore).
- 2) A 10mM aqueous HEPES-sodium hydroxide solution (pH 7.4), containing  $5 \times 10^{-3}\%$ (w/v) Tween 20, 0.20M sodium nitrate, and 10 $\mu$ M zinc nitrate was used as the running buffer. The sensor chip temperature was set at 25°C and the flow rate of the running buffer was set at 30 $\mu$ L/min. The running buffer was made to flow until the surface plasmon resonance value stabilized. This procedure was carried out on flow cell A.
- 3) As a biotin derivative of Phos-tag<sup>®</sup>, Phos-tag<sup>®</sup> BTL-105 was used and this was dissolved in the running buffer ( $5 \times 10^{-3}\%$ (w/v) Tween 20, 0.20M sodium nitrate, and 10 $\mu$ M zinc nitrate in 10mM aqueous HEPES-sodium hydroxide solution (pH 7.4)).
- 4) Binding onto flow cell A was carried out for 6 minutes at a ligand concentration of 1.0mM, a temperature of 25°C, and a flow rate of 30 $\mu$ L/min.
- 5) For the interaction of the Phos-tag<sup>®</sup> biotin and streptavidin, a surface plasmon resonance signal of 240RU was obtained for the maximum binding amount.
- 6) A sensor chip (Phos-tag<sup>®</sup> sensor chip), having a Phos-tag<sup>®</sup>-bound sensor, was prepared by the above procedures.

### Surface plasmon analysis of $\beta$ -casein

- 1)  $\beta$ -casein (pentaphosphorylated protein, Sigma Co., Ltd.) was used as the sample to be analyzed, and this sample was dissolved in the running buffer ( $5 \times 10^{-3}\%$ (w/v) Tween 20, 0.20M sodium nitrate, and 10 $\mu$ M zinc nitrate in 10mM aqueous HEPES-sodium hydroxide solution (pH 7.4)) (sample solution).
- 2) The surface plasmon resonance analysis was carried out in flow cell A at a sample concentration of 1.5 $\mu$ M, a temperature of 25°C, a flow rate of 30 $\mu$ L/min, a binding time of 15 minutes and a dissociation time of 10 minutes.
- 3) After measurement, the sensor chip was reactivated (removal of residual

bound substances) by the addition of a 400mM aqueous phosphate solution (pH 7.0) for 6 minutes, a 200mM aqueous EDTA (pH 8.0) solution for 6 minutes and the abovementioned running buffer for 5 minutes.

- 4) The measurement results are shown in FIG. 1.
- 5) The results of FIG. 1 shows that a high RU value is obtained for each flow cell and that the maximum binding amount for flow cell A is 2056.

### Interaction of the streptavidin sensor chip and $\beta$ -casein

- 1)  $\beta$ -casein (pentaphosphorylated protein, Sigma Co., Ltd.) was used as the sample to be analyzed, and this sample was dissolved in the running buffer (5X10<sup>-3</sup>%(w/v) Tween 20, 0.20M sodium nitrate, and 10 $\mu$ M zinc nitrate in 10mM aqueous HEPES-sodium hydroxide solution (pH 7.4)) (sample solution).
- 2) The surface plasmon resonance analysis was carried out in flow cell A of a new sensor chip at a sample concentration of 1.5 $\mu$ M, a temperature of 25°C, a flow rate of 30 $\mu$ L/min, a binding time of 15 minutes and a dissociation time of 10 minutes.
- 3) There was absolutely no interaction of  $\beta$ -casein with streptavidin.




FIG. 1. Surface plasmon analysis of phosphorylated  $\beta$ -casein using Phos-tag<sup>®</sup> Biotin (BTL-105) and Sensor Chip SA