(109 × 209mm Size)

FUJIFILM

Wako

Code No. 293-58301

SLP-HS Single Reagent Set

《Research Use only, Not for use
in diagnostic procedures》

Preface

The hemolymph of the silkworm *Bombyx mori* contains a self-defense mechanism termed the “prophenoloxidase cascade system”. Upon invasion of foreign material such as bacteria and fungi, the cascade system participates in melanin formation observed in the body fluid of insects to protect them from the invader’s attack. The system is triggered by peptidoglycan (PG) from bacteria and (1→3)- β -D-glucan (β -glucan) from fungi and yeast, and consequently, prophenoloxidase in the system is activated. It is also postulated that serial serine proteases are involved in the cascade system, but little has been elucidated on this point.

SLP (Silkworm Larvae Plasma) is prepared sterile and contains all of the factors participating in the activation of the cascade system. The system, containing L-3, 4-dihydroxyphenylalanine (DOPA) as a phenoloxidase substrate, is triggered by PG and/or β -glucan is widely distributed in the cell wall of a variety of fungi. Therefore, SLP enables one to quantify and/or detect contamination by bacterial and fungal material by measuring the development of melanin formation in the assay.

Features

1. High sensitive detection of bacteria and fungi by measuring PG and β -glucan
2. Accurate and sensitive quantification of PG and β -glucan by use of Toxinometer (Wako Pure Chemical Industries Ltd., Osaka, Japan)
3. Visual estimation of PG and β -glucan by color development without any special apparatus
4. This product is stable at 2-10°C for at least 24 months.

Assay Principle

The activation mechanism of SLP is considered as shown in Figure 1. Either PG or β -glucan binds to the respective recognition protein termed PGRP or β GRP. The prophenoloxidase cascade system is initiated by these reactions, and consequently, prophenoloxidase is converted to phenoloxidase. The phenoloxidase catalyzes oxidation of DOPA and is followed by the formation of melanin in the mixture.

One can visually judge a reaction by the change in color or the reaction mixture after a certain period of time due to the development of a block color.

When using Toxinometer, one can quantitate either PG or β -glucan by measuring the Activating time (T_a) of the reaction and the absorbance change caused by the color change from the start of the reaction to a predetermined level.


Fig 1. Activation mechanism of SLP system

Kit contents

1. SLP-HS(Silkworm Larvae Plasma High Sensitive)
for 0.2 mL \times 20 vials
Sensitivity : Either 10 pg/mL of PG or 1 pg/mL of β -glucan can be detected with Toxinometer at 30°C in 120 minutes.
2. SLP-Diluent 1.0 mL \times 20 vials
3. Standard (Digested Peptidoglycan from *S. aureus*)
0.5 mL \times 1 vial

Assay Procedures

<Procedure 1>

- 1) Dispense 0.1 mL of SLP-Diluent provided into the vial of SLP-HS, and then gently stir it to dissolve the SLP material completely to reconstitute the SLP reagent.

- 2) Add 0.1 mL of either dilution series of Standard provided or samples to the SLP reagent.
- 3) Mix the solutions by gently stirring the vial, set the vial to Toxinometer at 30 °C, and measure it.

〈**Procedure 2**〉

- 1) Add either the dilution series of Standard or samples to SLP-diluent vial to become the more than 2 volumes of dilution for added solution.
- 2) Remove 0.2 mL of diluted sample with SLP-Diluent to transfer it to a SLP-HS vial.
- 3) Dissolved SLP-HS completely with gently stirring the vial, set the vial to Toxinometer and measure it as described above.

※See the details in “Instruction of Toxinometer”.

※Visual estimation : Concentration of SLP activation material in the samples can be determined based on the endpoint and the dye development of the serial dilution solutions with the standard material.

Usage notes

1. All the vessels and apparatuses to be used should be free of SLP activating substance such as PG and β -glucan. Dry-heating at 250 °C for 2 hours is recommended.
2. The sensitive of SLP reagent is suppressed by samples with a high ionic strength.
3. This kit is for research use only and not for use in diagnostic procedures. For in vitro use only. Not for internal use in humans or animals.

Storage Keep at 2~10 °C

Expiration Indicate on the label (outer box)

Package 20 Tests

References

- 1) Ashida, M. and Yamazaki, H. I. : “Molting and Metamorphosis” ed. by Ohnishi, E. and Ishizaki, H., Japan Sci. Soc. Press, Tokyo/Springer-Verlag, Berlin, p.239 (1990)
- 2) Tsuchiya, M., Asahi, N., *et al.* : Detection of peptidoglycan and β -glucan with silkworm larvae plasma test. FEMS Immunol. Med. Microbiol., **15**, 129-134 (1996)

Related product

297-51501	SLP Reagent Set	for 3 mL
293-26551	Limulus Test Tube-S (12 × 75mm, Endotoxin Free)	10 vials × 10
293-28251	Aluminium Cap-S (14.7 × 18 mm, Endotoxin Free)	10 pieces × 10
298-35031	Bio Clean Tip Wako [®] 1000 II	100 tips
291-35021	Bio Clean Tip Wako [®] 200 II	100 tips
294-35011	Bio Clean Tip Wako [®] Extend S II	100 tips

FUJIFILM Wako Pure Chemical Corporation

1-2, Doshomachi 3-Chome, Chuo-Ku, Osaka 540-8605, Japan
Telephone : + 81-6-6203-3741
Facsimile : + 81-6-6201-5964
<http://www.wako-chem.co.jp>

FUJIFILM Wako Chemicals U.S.A. Corporation

1600 Bellwood Road
Richmond, VA 23237
U.S.A.
Telephone : + 1-804-271-7677
Facsimile : + 1-804-271-7791
<http://www.wakousa.com>

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