< For Research Use Only, Not for USE in Diagnostic Procedure >

Code No. 297-51501

< For Microorganism Detection >

SLP Reagent Set

INTENDED USE

Wako SLP Reagent Set is intended for detection of peptidoglycan and (1→3)- β -D-glucan.

SUMMARY AND GENERAL INFORMATION

The hemolymph of the silkworm *Bombyx mori* contains a self-defense mechanism termed the “prophenoloxidase cascade system”. Upon invasion of foreign material such as bacteria and fungi, the cascade system participates in melanin formation observed in the body fluids of insects to protect them from the invader’s attack. The system is triggered by peptidoglycan (PG) from bacteria and (1→3)- β -D-glucan (β -glucan) from fungi and yeast, and consequently, prophenoloxidase in the system is activated. It is also postulated that serial serine proteases are involved in the cascade system, but little has been elucidated on this point.

SLP (silkworm larvae plasma) is prepared sterile and contains all of the factors participating in the activation of the cascade system. The system, reconstituted with a solution containing L-3,4-dihydroxyphenylalanine (DOPA) as a phenoloxidase substrate, is triggered by PG and/or β -glucan. DOPA is oxidized by phenoloxidase, forming melanin pigment in the reaction mixture. PG is a component of most bacterial cell walls and β -glucan is widely distributed in the cell wall of a variety of fungi. Therefore, SLP enables one to quantify and/or detect contamination by bacterial and fungal material by measuring the development of melanin formation in the assay.

Features

1. Sensitive quantitation and detection of PG and β -glucan.
2. Sensitive and accurate quantitation of PG and β -glucan by use of a microplate reader or Toxinometer.
3. Visual detection of bacterial and fungal contamination without any special apparatus.
4. SLP in a lyophilized form is stable at 2~10°C for at least 18 months.

PRINCIPLE

The activation mechanism of SLP is considered as shown in Figure 1. Either PG or β -glucan binds to the respective recognition protein termed PGRP or β GRP. The prophenoloxidase cascade system is initiated by these reactions, and consequently, prophenoloxidase is converted to phenoloxidase. The phenoloxidase catalyzes oxidation of DOPA and is followed by the formation of melanin in the mixture.

One can visually judge reaction by the change in color of the reaction mixture after a certain period of time due to the development of a black color.

When using a microplate reader or Toxinometer, one can quantitate either PG or β -glucan by measuring the onset time (T_a) of the reaction and the absorbance change caused by the color change from the start of the reaction to a predetermined level.


Fig 1 (Activation mechanism of SLP system)

REAGENTS

All the reagents should be stored at 2~10°C.

1. SLP lyophilized 1 vial for 3 mL

For use with a microplate reader or for visual tests, approximately 50 assays can be made. For use with a Toxinometer, approximately 30 assays can be performed.

Sensitivity ; Endpoints determined by the visual test described below after 1 hour incubation at 30°C with PG and β -glucan are indicated in each lot.

2. Substrate, lyophilized 1 vial for 3 mL

Contains DOPA as a phenoloxidase substrate.

3. Diluent 1 vial 4 mL

Contains a "good's" buffer.

Preparation of reagents

- (A) Add 3 mL of Diluent to substrate. Gently stir the mixture to dissolve the substrate completely. Avoid getting the rubber cap wet.
- (B) Add the entire reconstituted substrate solution to the vial of lyophilized SLP. Gently stir to completely dissolve the lyophilized product. Again, avoid wetting the rubber cap. The resultant solution is termed the "SLP test solution".

PROCEDURE

1. Visual test using microplate

- (A) Make a 2-fold dilution series with the reference material of either PG or β -glucan and add 0.05 mL of each dilution and sample to the wells of the microplate.
- (B) Dispense 0.05 mL of SLP test solution to each well.
- (C) After gentle stirring, incubate the plate at 30°C for one hour.
- (D) Visually judge the formation of melanin in the reaction mixture : sample where melanin formation is observed are positive and others are negative.
- (E) The minimum concentration of the reference material at which melanin is formed is defined as the endpoint of the SLP test solution.
- (F) The determination of whether the sample is negative or positive depends on whether or not the sample contains SLP reactive materials at a concentration greater than that of the endpoint.

2. Procedure using microplate reader

- (A) 0.05 mL of samples and the dilution series of a reference material of either PG or β -glucan are dispensed into each well of the microplate.
- (B) Add 0.05 mL of SLP test solution to each well quickly. The use of a sterilized successive dispenser is recommended.
- (C) Read the absorbance of the microplate reader under the following conditions :

Assay absorbance :	650 nm
Assay mode :	Kinetic
Reaction time :	90 minutes
Assay temperature :	30°C
Auto mix :	Once
- (D) The standard curve is plotted with a regression of Log-Log : horizontal axis-concentration of SLP activation material ; vertical axis-onset time.
- (E) Using the standard curve, calculate the concentration of SLP reactive material from the measured onset time.

3. Procedure using Toxinometer

When using the Toxinometer for the SLP test assay, the procedure can be performed in the same manner as the endotoxin assay described in the operator's manual. The Toxinometer may exhibit a signal indicating "lack of transmittance light" in the activated sample reaction mixture. Ignore this alarm. The results are not affected by the alarm.

Notes in using SLP

1. Apparatuses should be sterilized by dry-heating at 250 °C for at least 2 hours and should be free of contamination by SLP reactive material.
2. After reconstitution, store SLP test solutions at 0~4 °C and use them as soon as possible.
3. Pay careful attention to avoid contamination of the apparatuses by SLP activation material during storage.
4. The sensitivity of the SLP test solution is largely influenced by the ionic strength.
5. This kit is for research use only and not for use in diagnostic procedures. For in vitro use only. Not for internal use in humans or animals.

FUJIFILM Wako Pure Chemical Corporation

1-2, Doshomachi 3-Chome, Chuo-Ku, Osaka 540-8605, Japan
Telephone : + 81-6-6203-3741
Facsimile : + 81-6-6201-5964
<http://www.wako-chem.co.jp>

FUJIFILM Wako Chemicals U.S.A. Corporation

1600 Bellwood Road
Richmond, VA 23237
U.S.A.
Telephone : + 1-804-271-7677
Facsimile : + 1-804-271-7791
<http://www.wakousa.com>

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