## Limulus HS-J Single Test *Wako*

25 tests

### Package insert

#### [Summary and General Information]

Endotoxin, lipopolysaccharide (LPS) is a component of the outer membrane of Gram-negative bacteria. Since endotoxin causes serious symptoms such as fever and shock in the human body, the bacterial endotoxin test for pharmaceuticals, especially parenteral drugs, is important.

In 1956, F. B. Bang reported clotting of the horseshoe crab body fluid due to Gram negative bacteria<sup>1)</sup>. In 1964, J. Levin and F. B. Bang discovered that horseshoe crab amebocytes (Limulus amebocyte lysate, or LAL) contained all the components needed to clot due to endotoxin<sup>2)</sup>. Subsequently, bacterial endotoxin detection by the Limulus test has been widely used. This method known as the Bacterial Endotoxin Test was recognized in 1980 by United States Pharmacopoeia and, in 1988 by Japanese Pharmacopoeia<sup>3) 4)</sup>.

This kit contains LAL reagent and Control Standard Endotoxin (CSE) and can rapidly assay and detect endotoxin with high sensitivity. In addition, the labelled gel-clot sensitivity was determined using the Japanese Pharmacopoeia Reference Standard Endotoxin (JP-RSE). CSE contains 500 ng LPS (Lipopolysaccharide, LPS) which has been extracted and purified from *E. coli* UKT-B by phenol method<sup>5) 6)</sup>, and combined with mannitol and glycine as additives, and is tested in comparison to the Japanese Pharmacopoeia Reference Standard Endotoxin and the potency is printed on the label as a reference.

#### [Features]

1. Highly sensitive for detection of endotoxin.
2. Gel forms strongly and easy to determine.
3. Single test vials mean there is no reagent wasted.
4. Store at 2~10°C.

#### [Principle]

In the gelation of LAL by endotoxin, serine proteases are activated in the cascade reactions. In the last reaction, the gel-forming protein precursor (coagulogen) is hydrolyzed to coagulin, which makes an insoluble gel (Fig. 1).

The gel-clot technique is that LAL reagent is mixed with a sample. The mixture is incubated at 37°C for one hour. It is assayed by turning the tube upside down to determine if a gel has formed.


Fig. 1 the cascade system of Limulus Amebocyte Lysate

**[Kit Contents]**

1. LAL Reagent, lyophilized 25 tubes containing lyophilized extract of *Limulus polyphemus* amoebocytes  
..... 25 vials (for 0.20 mL)

Gel-clot sensitivity was determined using the Japanese Pharmacopoeia Reference Standard Endotoxin

Store at 2 to 10°C

2. Control Standard Endotoxin, lyophilized 1 vial × 500 ng (as purified LPS) Endotoxin purified from *E. coli* UKT-B, containing mannitol and glycine as additives.

Store at 2 to 10°C

**[Usage]**

I. Materials required but not supplied

1. Pipettes (for 0.2 mL, 0.5 mL, 2.0 mL and 5.0 mL)
2. Test tubes with aluminum caps for dilution of standard solutions
3. Block heater or an incubator able to maintain 37 ± 1°C.
4. Endotoxin-free water

NOTE : A test tube and the pipette should be dehydrogenated by 250°C and more than 30 min. Water should be Endotoxin-free water. The plastic products such as tips for micropipets should be confirmed that they are no endotoxin contamination and no interference for the measurement.



II. Procedure

1. Preparation of Control Standard Endotoxin (CSE) Solution  
Remove the aluminum cap on the vial. The inside of the vial is a vacuum. Release the vacuum by pulling slowly up on the rubber stopper and remove the stopper from the vial. The EU value (EU/vial) is printed on the label. Determine the volume of water to be added into the vial to make a 1,000 EU/mL solution. Add the determined volume of endotoxin-free water into the vial. Put the rubber stopper on the vial. After turning the vial upside down several times, agitate the vial vigorously with a Vortex mixer for approximately 2 minutes. The Control Standard Endotoxin Solution can be stored at 2~10°C for 1 month. When the stored solution is used, agitate the solution vigorously with Vortex mixer for more than 1 minute before use.
2. Dilution of standard endotoxin solutions, the pH adjustment and dilution of samples  
A dilute standard solution should be agitated with a Vortex mixer for more than 30 seconds before use in the next dilution procedure. One dilution should not be more than 10-fold.  
When the pH of the LAL reagent solution reconstituted with 0.20 mL of sample solution is out of the range between 6.0 and 8.0, adjust the pH range of sample solution with a dilute sodium hydroxide solution or dilute hydrochloric acid.
3. Dispensing of samples  
Take LAL reagent tubes out of the package and put the tubes in a tube rack. When powder of LAL is attached to the rubber stopper, knock on several times of bottoms of the vial. A little LAL attaches to a rubber stopper and does not influence the result of the examination. Each rubber stopper is loosened by releasing the vacuum and slowly pulling up the rubber stopper slightly. (Fig. 2)

Dispense 0.20 mL of a sample into each LAL reagent tube and put the rubber stopper on the tube. Mix gently to dissolve the LAL reagent without making bubbles. Incubate the reaction mixture at  $37 \pm 1^\circ\text{C}$  for  $60 \pm 2$  min without vibration. (Fig. 3)

4. Judgment

After 60 minute incubation, slowly turn each tube upside down to determine if a gel is formed. A positive result is defined as the formation of a firm gel capable of maintaining its integrity. A negative result is defined as the absence of gel or the formation of a viscous mass, which does not hold. (Fig. 4)

III. Interfering Testing

The inhibition or enhancement of the LAL test should be determined for each sample formulation before the LAL test is used to assess the endotoxin content of any sample. When the pH of the LAL reagent solution reconstituted with 0.20 mL of sample solution is out of the range between 6.0 and 8.0, adjust the pH range of sample solution with a dilute sodium hydroxide solution or dilute hydrochloric acid.

(Method)

1. Prepare a series of two-fold dilutions of endotoxin with endotoxin-free water to give concentrations of  $0.0078 \sim 0.0625$  EU/mL. Add  $50 \mu\text{L}$  of  $6.25$  EU/mL endotoxin solution to  $4.95$  mL of sample solution, which makes a final endotoxin concentration of  $0.0625$  EU/mL. Prepare another series of two-fold dilutions of the  $0.0078$  EU/mL-spiked sample solution diluted with the sample solution to give final endotoxin concentrations of  $0.0078 \sim 0.0625$  EU/mL.
2. Perform the test on the above 2 series of dilutions, an unspiked sample solution and a negative control (endotoxin-free water). Please follow "Inhibition and Enhancement Testing" of the station "bacterial endotoxins test" in the Japanese Pharmacopoeia about the number of the repetition.
3. The unspiked sample solutions and the negative controls must be negative. The geometric mean endpoint concentration of the first series of dilutions (diluted with endotoxin-free water) must be greater than or equal to  $0.0156$  EU/mL and less than or equal to  $0.0625$  EU/mL.  
If the negative controls are negative and the unspiked sample solutions are positive, remove endogenous endotoxin by ultrafiltration/appropriate way or dilute the sample solution. And then retest.

If all the conditions in the above 3. are met and the geometric mean endpoint concentration of the second series of dilutions (diluted with sample solution) is greater than or equal to 0.0156 EU/mL and less than or equal to 0.0625 EU/mL, it is confirmed that the sample concentration in the sample solution does not interfere with the LAL test.

#### IV. Preliminary Test for Dilution Factor

Prepare a series of dilutions (determine an appropriate dilution factor according to the degree of interference) of the sample solution with endotoxin-free water. Add endotoxin solution to each of the dilutions making the final endotoxin concentration  $2\lambda$  (0.0625 EU/mL). Perform the test on 2 series of dilutions (a series of unspiked sample solutions and a series of  $2\lambda$  endotoxin-spiked sample solutions). If a dilution of the spiked sample solution is positive and the same dilution of the unspiked sample solution is negative, it is confirmed that no interference is observed the sample concentration of that dilution.

#### V. Routine Testing (Limit Test)

Perform the test on the following specimens at least in duplicate (n=2).

- Negative control (Endotoxin free water)
- A dilution of endotoxin ( $2\lambda$ ) (Positive Control : PC)
- Sample solution
- Positive product control ( $2\lambda$  endotoxin-spiked sample solution) (PPC)

If all the following conditions are met, the endotoxin concentration of the sample solution is passed (The sample is less than 0.03 EU/mL).

- Negative control is negative.
- PC and PPC are positive
- Sample is negative

If the sample solution is positive, perform the test on a series of 2-fold dilutions of sample solution diluted with endotoxin-free water.

The endotoxin concentration in the sample solution is determined by :

$[\lambda \text{ (EU/mL)}] \times [\text{Maximum dilution factor in the series of dilutions of positive results}]$

#### **[Precautions]**

1. Be very careful not to cause endotoxin contamination. Use endotoxin free equipment and water.
2. Do not use the reagent with a large quantity of insoluble material or when it changes color.
3. For research used only. For in vitro use only. Not for use in diagnostic procedures. Not for internal use in human or animals.
4. Exercise caution when handling LAL, because its toxicity is not known.

**[Procedure to Determine Potency of CSE (Reference)]**

1. Prepare a series of dilutions of Japanese Pharmacopoeia Reference Standard Endotoxin (JP-RSE) and 4 vials of CSE according to Table 1 and 2. Perform the gel-clot technique test on the dilutions of 0.0078~0.0625 EU/mL (JP-RSE) and 0.00078~0.00625 ng/mL in tetra-plicate.

Calculate endpoint geometric mean of JP-RSE (refer to Table3) and CSE (refer to Table4). Calculate Endotoxin Unit per CSE weight (EU/ng). CSE potency (EU/vial) is multiplied 500 (CSE contains 500 ng).

The potency which calculated with the lot of LAL reagent and CSE is printed on the label as a reference.

Note : If a LAL lot is changed, the CSE potency test must be repeated because EU/ng of the CSE may be changed.

Table 1. Endotoxin dilution series of JP-RSE

JP-RSE soln. (mL)	H <sub>2</sub> O (mL)	final JP-RSE conc. (EU/mL)
0.5 (10000 EU/mL)	4.5	1000
0.5 (1000 EU/mL)	4.5	100
0.5 (100 EU/mL)	4.5	10
0.5 (10 EU/mL)	4.5	1.0
1.0 (1.0 EU/mL)	3.0	0.25
1.0 (0.25 EU/mL)	3.0	0.0625
2.0 (0.0625 EU/mL)	2.0	0.0313
2.0 (0.0313 EU/mL)	2.0	0.0156
2.0 (0.0156 EU/mL)	2.0	0.0078

Table 2. Endotoxin dilution series of CSE

CSE soln. (mL)	H <sub>2</sub> O (mL)	final CSE conc. (ng/mL)
0.5 (100 ng/mL)	4.5	10
0.5 (10 ng/mL)	4.5	1.0
0.5 (1.0 ng/mL)	4.5	0.1
1.0 (0.1 ng/mL)	3.0	0.025
1.0 (0.025 ng/mL)	3.0	0.00625
2.0 (0.00625 ng/mL)	2.0	0.00313
2.0 (0.00313 ng/mL)	2.0	0.00156
2.0 (0.00156 ng/mL)	2.0	0.00078

Table 3 Example of JP-RSE endpoint

Replication	NC	0.0625	0.0313	0.0156	0.0078	endpoint(EU/mL)
1	-	+	+	+	-	0.0156
2	-	+	+	-	-	0.0313
3	-	+	+	-	-	0.0313
4	-	+	+	-	-	0.0313

Table 4 Example of CSE endpoint (Vial 1)

Replication	0.00625	0.00313	0.00156	0.00078	endpoint(ng/mL)
1	+	+	+	-	0.00156
2	+	+	-	-	0.00313
3	+	+	-	-	0.00313
4	+	+	-	-	0.00313

CSE (Vial 2)

1	+	+	-	-	0.00313
2	+	+	-	-	0.00313
3	+	+	-	-	0.00313
4	+	+	-	-	0.00313

CSE (Vial 3)

1	+	+	-	-	0.00313
2	+	+	-	-	0.00313
3	+	+	-	-	0.00313
4	+	-	-	-	0.00625

CSE (Vial 4)

1	+	+	-	-	0.00313
2	+	+	-	-	0.00313
3	+	+	+	-	0.00156
4	+	+	-	-	0.00313

Geometric Mean (GM) of JP-RSE

$$(\log(0.0156) + 3\log(0.0313))/4 = -1.580$$

$$\text{antilog}(-1.580) = 0.0263 \text{ EU/mL}$$

GM of CSE

$$(2\log(0.00156) + 13\log(0.00313) + \log(0.00625))/16 = -2.523$$

$$\text{antilog}(-2.523) = 0.00300 \text{ ng/mL}$$

$$\text{CSE (EU/ng)} = 0.0263/0.00300 = 8.77 \text{ EU/ng} = 8.8 \text{ EU/ng}$$

$$\text{Therefore CSE potency is } 8.8 \text{ (EU/ng)} \times 500 \text{ (ng/vial)}$$

$$= 4400 \text{ EU/vial}$$

### [References]

1. Bang, F. B. : *Bull. Johns Hopkins Hosp.*, **98**, 325 (1956).
2. Levin, J. and Bang, F. B. : *Bull. Johns Hopkins Hosp.*, **115**, 265 (1964).
3. *The United States Pharmacopoeia 22th, The National Formulary 17th*, Supplement 8, p.3349, U.S. Pharmacopoeial Convention Inc. (1990).
4. THE JAPANESE PHARMACOPOEIA TWELVE EDITION Handbook, B-57 (1991)
5. Westphal, O., Lüderitz, O., Eichenberger, E. and Keiderling, W. : *Z. Naturforsch.*, **7b**, 536 (1952).
6. Akama, K., Kuratsuka, K., Homma, R., Kanoh, S., Niwa, M., Iwanaga, S. and Nakahara, C. : *Bacterial Endotoxin*, Eds. Homma, J.Y. *et al.*, p.395, Verlag Chemie (1984).

---

## FUJIFILM Wako Pure Chemical Corporation

1-2, Doshomachi 3-Chome, Chuo-Ku, Osaka 540-8605, Japan  
Telephone : + 81-6-6203-3741  
Facsimile : + 81-6-6201-5964  
<http://www.wako-chem.co.jp>

### FUJIFILM Wako Chemicals U.S.A. Corporation

1600 Bellwood Road  
Richmond, VA 23237  
U.S.A.  
Telephone : + 1-804-271-7677  
Facsimile : + 1-804-271-7791  
<http://www.wakousa.com>

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Laboratories, Inc.