

Wako**BioProducts****Anti Serum for Immunochemistry**

<For Research Use Only>

These polyclonal antibodies, produced by rabbit or sheep using purified antigen immunization, and are respectively specific reactive to corresponding antigens in a sera and tissues of human.

The titers are determined by quantitative precipitin formation and are indicated by the amount of antibodies in 1 mL of anti serum. The lot variation of each serum is within 10 %. Each product contains 0.1 % NaN₃ as a preservative.

Code No.	Description		Titer	Package
014-10031	Anti Human C1q	Rabbit	-	1 mL
019-10221	Anti Human C3	Sheep	1,500	1 mL
012-10071	Anti Human C5	Rabbit	-	1 mL
019-10081	Anti Human C9	Rabbit	-	1 mL
014-10151	Anti C-Reactive Protein	Sheep	4,000	2 mL
016-10111	Anti Human Factor H	Rabbit	-	1 mL
013-10121	Anti Human Factor I	Rabbit	-	1 mL
013-10241	Anti Human Fibronectin	Rabbit	3,100	1 mL
018-10171	Anti Human -Lipoprotein	Sheep	-	2 mL
016-10091	Anti Human Properdin	Rabbit	-	1 mL

Storage: Keep at 2 – 10

020108IBD

Wako Pure Chemical Industries, Ltd.
1-2, Doshomachi 3-Chome
Chuo-Ku, Osaka 540-8605, Japan
Telephone: 81-6-6203-3741
Facsimile: 81-6-6201-5964
Online Cat.: <http://search.wako-chem.com>

Wako Chemicals USA, Inc.
1600 Bellwood Road
Richmond, VA 23237, U.S.A.
Toll-Free (U.S. only): +1-877-714-1920
Telephone: 1-804-714-1920
Facsimile: 1-804-271-7791
<http://www.wakousa.com>

Wako Chemicals GmbH
Nissanstraße 2, D-41468
Neuss, Germany
Telephone: 49-2131-311-0
Facsimile: 49-2131-311100
E-mail: biochem@wako.jis.de