



Antibody

---

## Anti GM3 (Neu Ac), Monoclonal Antibody

---

<For Research Use Only>

**Wako Catalog No. 303-06541 1.0 mL**

**309-06543 0.5 mL**

Storage Condition: Keep at 2 ~ 10°C

manufactured by Nihon Biotest Institute (Japan).

Monoclonal antibody to GM3 (Neu Ac) is the product of a hybridoma derived from the fusion of mouse splenocytes from BALB/c mice immunized with highly purified GM3 (Neu Ac) and myeloma cells (P3U1) by polyethylene glycol method.

This antibody is 10 mg/mL solution dissolved in 10 mM phosphate buffer (pH 7.2) containing 0.15 M NaCl.

**Preparation:** Purified from ascites fluid of BALB/c mice by ammonium sulfate fractionation and DEAE-cellulose column chromatography.

**Clone No.:** M2590

**Isotype:** IgM

---

Wako Pure Chemical Industries, Ltd.  
1-2, Doshomachi 3-Chome  
Chuo-Ku, Osaka 540-8605, Japan  
Telephone: 81-6-6203-3741  
Facsimile: 81-6-6201-5964  
Online Cat.: <http://search.wako-chem.com>

Wako Chemicals USA, Inc.  
1600 Bellwood Road  
Richmond, VA 23237, U.S.A.  
Toll-Free (U.S. only): +1-877-714-1920  
Telephone: 1-804-714-1920  
Facsimile: 1-804-271-7791  
<http://www.wakousa.com>

Wako Chemicals GmbH  
Nissanstraße 2, D-41468  
Neuss, Germany  
Telephone: 49-2131-311-0  
Facsimile: 49-2131-311100  
<http://www.wakochemicals.de>