

品名	ページ
A	
Acrp30, Globular Domain, Mouse, recombinant (gAcrp30)	76
Activin A, Human, recombinant, Animal-derived-free	58
AITRL, Human, recombinant	25
Amphiregulin, Human, recombinant	37
Apolipoprotein A-1, Human, recombinant (ApoA-1)	78
Apolipoprotein E2, Human, recombinant (ApoE2)	78
Apolipoprotein E3, Human, recombinant (ApoE3)	78
Apolipoprotein E4, Human, recombinant (ApoE4)	78
Apo-SAA1, Human, recombinant	79
APRIL, Human, Insect Cells recombinant	25
Artemin, Human, recombinant	48
B	
BAFF, Human, recombinant	25
4-1BBL soluble, Human, recombinant (CD157)	25
4-1BB Receptor soluble, Human, recombinant (4-1BBR/CD137)	29
Betacellulin, Human, recombinant	35
Betacellulin, Mouse, recombinant	35
Betacellulin, Rat, recombinant	35
BLC/BCA-1, Mouse, recombinant (CXCL13)	69
Bone Morphogenetic Protein 2, Human, recombinant (BMP-2)	58
Bone Morphogenetic Protein 3, Human, recombinant (BMP-3)	58
Bone Morphogenetic Protein 4, Human, recombinant (BMP-4)	58
Bone Morphogenetic Protein 4(truncated), Human, recombinant (BMP-4(truncated))	59
Bone Morphogenetic Protein 6, Human, recombinant (BMP-6)	59
Bone Morphogenetic Protein 7, Human, recombinant (BMP-7)	59
Bone Morphogenetic Protein 13, Human, recombinant (BMP-13)	59
Brain Derived Neurotrophic Factor, Human, recombinant (BDNF)	46
Brain Derived Neurotrophic Factor, Human, recombinant, Animal-derived-free (BDNF)	46
C	
C10, Mouse, recombinant (CCL6/MRP-1)	71
Cardiotrophin-1, Human, recombinant	79
Cardiotrophin-1, Mouse, recombinant	79
CD30 Ligand soluble, Human, recombinant (sCD30L/CD153)	26
CD40 Ligand soluble, Human, recombinant (sCD40L/TNFSF5)	26
Ciliary Neurotrophic Factor, Human, recombinant (CNTF)	49
CTACK, Human, recombinant (CCL27)	75
CTACK, Mouse, recombinant (CCL27)	75
CTGF, Human, recombinant (CCN2)	50
CTGFL, Human, recombinant (WISP-2/CCN5)	50
D	
DLL-1 soluble, Human, recombinant	64
E	
Epidermal Growth Factor, Human, recombinant, Animal-derived-free (EGF)	35
Epidermal Growth Factor, Mouse, recombinant (EGF)	36
Epidermal Growth Factor, Rat, recombinant (EGF)	36
EGF Receptor soluble, Human, recombinant	36

品名	ページ
EMAP-II, Human, recombinant (AIMP/SCYE1)	79
ENA-78, Human, recombinant (CXCL5)	67
Eotaxin-2, Mouse, recombinant (CCL24)	74
Eotaxin-3, Human, recombinant (CCL26)	75
Epigen, Human, recombinant	36
Epiregulin, Human, recombinant	36
E-Selectin, Human, recombinant (CD62E)	79
Exodus-2/SLC, Human, recombinant (CCL21)	73
Exodus-2/SLC, Mouse, recombinant (CCL21)	74
F	
Fas Ligand soluble, Human, recombinant (FASL/CD95L)	26
Fibroblast Growth Factor (acidic), Human, recombinant (acidic FGF/FGF1)	37
Fibroblast Growth Factor (acidic), Human, recombinant, Animal-derived-free (acidic FGF/FGF1)	38
Fibroblast Growth Factor (basic), Mouse, recombinant (basic FGF/FGF2)	38
Fibroblast Growth Factor (basic), Mouse, recombinant, Animal-derived-free (basic FGF/FGF2)	39
Fibroblast Growth Factor (basic), Human, recombinant (basic FGF/FGF2)	38
Fibroblast Growth Factor (basic), Human, recombinant (basic FGF/FGF2)	38
Fibroblast Growth Factor (basic), Human, recombinant, Animal-derived-free (basic FGF/FGF2)	38
Fibroblast Growth Factor 4, Human, recombinant (FGF4)	39
Fibroblast Growth Factor 4, Human, recombinant, Animal-derived-free (FGF4)	39
Fibroblast Growth Factor 5, Human, recombinant (FGF5)	39
Fibroblast Growth Factor 6, Human, recombinant (FGF6)	40
Fibroblast Growth Factor 8, Human, recombinant (FGF8)	40
Fibroblast Growth Factor 9, Human, recombinant (FGF9)	41
Fibroblast Growth Factor 9, Mouse, recombinant (FGF9)	41
Fibroblast Growth Factor 10, Human, recombinant (FGF10)	41
Fibroblast Growth Factor 10, Human, recombinant, Animal-derived-free (FGF10)	41
Fibroblast Growth Factor 16, Human, recombinant (FGF16)	41
Fibroblast Growth Factor 17, Human, recombinant (FGF17)	42
Fibroblast Growth Factor 19, Human, recombinant (FGF19)	42
Fibroblast Growth Factor 20, Human, recombinant (FGF20)	42
Fibroblast Growth Factor 21, Human, recombinant (FGF21)	42
Fibroblast Growth Factor 23, Human, recombinant (FGF23)	43
Flt3 Ligand, Human, recombinant (Flt3L)	20
Flt3 Ligand, Human, recombinant, Animal-derived-free (Flt3L)	21
Flt3 Ligand, Mouse, recombinant (Flt3L)	21
Follistatin, Human, recombinant	62
Fractalkine (CX3CL1)	75
sFRP-1, Human, recombinant	64
G	
GDF-2, Human, recombinant	59
GDF-3, Human, recombinant	60
GDF-7, Human, recombinant	60
GDF-11, Human, recombinant	60
Glial Cell Line-derived Neurotrophic Factor, Human, recombinant (GDNF)	48
Glial Cell Line-derived Neurotrophic Factor, Mouse, recombinant (GDNF)	48
Glial Cell Line-derived Neurotrophic Factor, Rat, recombinant (GDNF)	49
GMF- β , Human, recombinant	49
Granulocyte Colony-Stimulating Factor, Human, recombinant (G-CSF)	18

品名	ページ
Granulocyte Colony-Stimulating Factor, Human, recombinant, Animal-derived-free (G-CSF)	18
Granulocyte Colony-Stimulating Factor, Mouse, recombinant (G-CSF)	18
Granulocyte-Macrophage Colony-Stimulating Factor, Human, recombinant (GM-CSF)	18
Granulocyte-Macrophage Colony-Stimulating Factor, Human, recombinant, Animal-derived-free (GM-CSF)	19
Granulocyte-Macrophage Colony-Stimulating Factor, Mouse, recombinant (GM-CSF)	19
Granulocyte-Macrophage Colony-Stimulating Factor, Mouse, recombinant, Animal-derived-free (GM-CSF)	19
Granulocyte-Macrophage Colony Stimulating Factor, Rat, recombinant (GM-CSF)	19
GRO, Rat, recombinant (CXCL1/KC)	66
GRO - α , Human, recombinant (CXCL1/MGSA)	66
GRO - β , Rat, recombinant (CXCL2/MIP-2)	66
GRO - γ , Human, recombinant (CXCL3)	66
H	
Hepatocyte Growth Factor, Human, Insect Cells recombinant (HGF)	51
Hepatocyte Growth Factor, Mouse, recombinant (expressed in Insect Cells) (HGF)	51
Hepatocyte Growth Factor, Rat, recombinant, Solution (HGF)	51
Heregulin - β -1, Human, Animal-derived-free (HRG - β 1)	52
HVEM-Fc, Human, recombinant(expressed in Insect Cells)	29
I	
I-309, Human, recombinant (CCL1)	69
ICAM-1, Human, recombinant	80
Insulin-like Growth Factor Binding Protein 1, Human, recombinant (IGF-BP1)	53
Insulin-like Growth Factor Binding Protein 2, Human, recombinant (expressed in Insect Cells) (IGF-BP2)	53
Insulin-like Growth Factor Binding Protein 3, Human, recombinant (IGF-BP3)	53
Insulin-like Growth Factor Binding Protein 4, Human, recombinant (expressed in Insect Cells) (IGF-BP4)	54
Insulin-like Growth Factor Binding Protein 5, Human, recombinant (IGF-BP5)	54
Insulin-like Growth Factor Binding Protein 6, Human, recombinant (expressed in Insect Cells) (IGF-BP6)	54
Insulin-like Growth Factor Binding Protein 7, Human, recombinant (IGF-BP7)	54
Insulin-like Growth Factor-I, Human, recombinant (IGF-I)	52
Insulin-like Growth factor-I, Human, recombinant, Animal-derived-free (IGF-I)	52
Insulin-like Growth Factor-I, Mouse, recombinant (IGF-I)	52
Insulin-like Growth Factor-II, Human, recombinant (IGF-II)	53
Insulin-like Growth Factor-II, Human, recombinant, Animal-derived-free (IGF-II)	53
Interferon - β , Human, recombinant (IFN - β)	15
Interferon - γ , Human, recombinant (IFN - γ)	15
Interferon - γ , Human, recombinant, Animal-derived-free (IFN - γ)	15
Interferon - γ , Mouse, recombinant (IFN - γ)	15
Interferon - γ , Rat, recombinant (IFN - γ)	16
Interferon - λ 1, Human, recombinant (IFN - λ 1)	16
Interferon - λ 2, Human, recombinant (IFN - λ 2)	16
Interferon - λ 2, Mouse, recombinant (IFN - λ 2)	16
Interferon - ω , Human, recombinant (IFN - ω)	16
Interleukin -1 α , Human, recombinant (IL -1 α)	2
Interleukin -1 α , Mouse, recombinant (IL -1 α)	2
Interleukin -1 α , Rat, recombinant (IL -1 α)	2
Interleukin -1 β , Human, recombinant (IL -1 β)	2
Interleukin -1 β , Human, recombinant, Animal-derived-free (IL-1 β)	3
Interleukin -1 β , Mouse, recombinant (IL -1 β)	3
Interleukin -1 β , Rat, recombinant (IL -1 β)	3
Interleukin-1Receptor Antagonist, Human, recombinant (IL-1RA)	12

品名	ページ
Interleukin-2, Human, recombinant (IL-2)	3
Interleukin-2, Human, recombinant, Animal-derived-free (IL-2)	3
Interleukin-2, mouse, recombinant (IL-2)	3
Interleukin-2, Rat, recombinant (IL-2)	4
Interleukin-2, Receptor α soluble, Human, Insect Cells (sIL-2R α)	13
Interleukin-3, Human, recombinant (IL-3)	4
Interleukin-3, Human, recombinant, Animal-derived-free (IL-3)	4
Interleukin-3, Mouse, recombinant (IL-3)	4
Interleukin-3, Mouse, recombinant, Animal-derived-free (IL-3)	4
Interleukin-3 β , Rat, recombinant (IL-3 β)	4
Interleukin-4, Human, recombinant (IL-4)	5
Interleukin-4, Human, recombinant, Animal-derived-free (IL-4)	5
Interleukin-4, Mouse, recombinant (IL-4)	5
Interleukin-4, Mouse, recombinant, Animal-derived-free (IL-4)	5
Interleukin-4, Rat, recombinant (IL-4)	5
Interleukin-4, Receptor α soluble, Human, recombinant (sIL-4R α)	13
Interleukin-5, Human, recombinant (IL-5)	6
Interleukin-6, Human, recombinant (IL-6)	6
Interleukin-6, Human, recombinant, Animal-derived-free (IL-6)	6
Interleukin-6, Mouse, recombinant (IL-6)	6
Interleukin-6, Rat, recombinant (IL-6)	7
Interleukin-6 Receptor α soluble, Human, recombinant (sIL-6R α)	13
Interleukin-7, Human, recombinant (IL-7)	7
Interleukin-7, Human, recombinant, Animal-derived-free (IL-7)	7
Interleukin-7, Mouse, recombinant (IL-7)	7
Interleukin-8 (endothelial cell-derived), Human, recombinant (IL-8)	7
Interleukin-8 (monocyte-derived), Human, recombinant (IL-8)	8
Interleukin-9, Human, recombinant (IL-9)	8
Interleukin-10, Human, recombinant (IL-10)	8
Interleukin-10, Mouse, recombinant (IL-10)	8
Interleukin-10, Rat, recombinant (IL-10)	8
Interleukin-11, Human, recombinant (IL-11)	9
Interleukin-12, Human, recombinant (IL-12)	9
Interleukin-12, Mouse, recombinant (IL-12)	9
Interleukin-12 p40, Human, recombinant (IL-12 p40)	9
Interleukin-12 p40, Mouse, recombinant (IL-12 p40)	9
Interleukin-13, Human, recombinant (IL-13)	10
Interleukin-13, Mouse, recombinant (IL-13)	10
Interleukin-13 Variant, Human, recombinant (IL-13 Variant)	10
Interleukin-16, Human, recombinant, Animal-derived-free (IL-16)	10
Interleukin-17, Human, recombinant (IL-17)	10
Interleukin-17F, Human, recombinant (IL-17F)	11
Interleukin-19, Human, recombinant (IL-19)	11
Interleukin-20, Human, recombinant (IL-20)	11
Interleukin-22, Human, recombinant (IL-22)	11
Interleukin-23, Human, recombinant (expressed in Insect Cells), Solution (IL-23)	11
Interleukin-23, Mouse, recombinant (expressed in Insect Cells), Solution (IL-23)	12
Interleukin-33, Human, recombinant (IL-33)	12
Interleukin-31, Mouse, recombinant (IL-31)	12
IP-10, Human, recombinant (CXCL10)	67
I-TAC, Human, recombinant (CXCL11)	67

品名	ページ
K	
Keratinocyte Growth Factor, Human, recombinant (KGF/FGF7)	40
Keratinocyte Growth Factor, Human, recombinant, Animal-derived-free (KGF/FGF7)	40
L	
LD78 β , Human, recombinant (CCL3L1/MIP-1)	69
Leptin, Mouse, recombinant (OB)	76
M	
Macrophage Colony Stimulating Factor, Human, recombinant (M-CSF)	19
Macrophage Colony-Stimulating Factor, Human, recombinant, Animal-derived-free (M-CSF)	20
Macrophage Colony Stimulating Factor, Mouse, recombinant (M-CSF)	20
Macrophage Colony-Stimulating Factor, Mouse, recombinant, Animal-derived-free (M-CSF)	20
Macrophage Colony Stimulating Factor, Rat, recombinant (M-CSF)	20
Maspin, Human, recombinant	80
MCP-1, Human, recombinant (CCL2)	70
MCP-1, Human, recombinant, Animal-derived-free (CCL2)	70
MCP-1, Mouse, recombinant (CCL2/JE)	70
MCP-1, Rat, recombinant (CCL2)	70
MCP-2, Human, recombinant (CCL8)	72
MCP-2, Mouse, recombinant (CCL8)	72
MCP-3, Human, recombinant (CCL7)	72
MCP-4, Human, recombinant (CCL13)	72
MDC, Human, recombinant (CCL22)	74
MDC, Mouse, recombinant (CCL22)	74
MEC, Human, recombinant (CCL28)	75
Midkine, Human, recombinant	49
MIG, Human, recombinant (CXCL9)	69
MIP-1 α , Human, recombinant (CCL3)	70
MIP-1 α , Mouse, recombinant (CCL3)	71
MIP-1 α , Rat, recombinant (CCL3)	71
MIP-1 β , Human, recombinant (CCL4)	71
MIP-2, Mouse, recombinant (CXCL2)	69
MIP-3 α /LARC, Human, recombinant (CCL20)	73
MIP-3 β /ELC, Human, recombinant (CCL19)	73
MIP-4/PARC, Human, recombinant (CCL18)	73
Myostatin, Human, recombinant (GDF-8)	60
N	
Nanog, Human, recombinant	64
NAP-2, Human, recombinant (CXCL7)	67
Nerve Growth Factor - β , Human, recombinant (NGF- β)	46
Nerve Growth Factor - β , Human, recombinant, Animal-derived-free (NGF- β)	47
Nerve Growth Factor - β , Mouse, recombinant (NGF- β)	47
Neuroserpin, Human, recombinant	50
Neurotrophin-3, Human, recombinant (NT-3)	47
Neurotrophin-3, Human, recombinant, Animal-derived-free (NT-3)	47
Neurotrophin-4, Human, recombinant (NT-4)	48
Noggin (Dimer), Human, recombinant (expressed in HEK293 Cells)	62
Noggin, Mouse, recombinant	62

	ページ
O	
Oncostatin M, Human, recombinant (OSM)	54
Oncostatin M, Rat, recombinant (OSM)	55
Osteonectin, Human, recombinant (SPARC/BM-40)	55
Osteopontin, Human, recombinant (OPN)	80
Osteoprotegerin(22-202), Human, recombinant (OPG)	29
OX40 Ligand, Human, recombinant (expressed in Insect Cells) (OX40L/CD252)	26
P	
Parathyroid Hormone Related Protein, Human, recombinant (PTHrP)	80
PDGF-AA, Human, recombinant	55
PDGF-AA, Human, recombinant, Animal-derived-free	55
PDGF-AA, Mouse, recombinant	55
PDGF-AB, Human, recombinant	56
PDGF-BB, Human, recombinant	56
PDGF-BB, Human, recombinant, Animal-derived-free	56
PDGF-BB, Mouse, recombinant	56
PECAM-1, Human, recombinant	81
PEDF, Human, recombinant	81
Persephin, Human, recombinant	49
PF-4, Human, recombinant (CXCL4)	66
Placenta Growth Factor-1, Human, recombinant, Animal-derived-free (PLGF-1)	43
Placenta Growth Factor-3, Human, recombinant (PLGF-3)	43
Pleiotrophin, Human, recombinant (PTN)	50
Prolactin, Human, recombinant	81
R	
RANK Ligand soluble, Human, recombinant (sRANKL)	27
RANK Ligand soluble, Mouse, recombinant (sRANKL)	27
RANK Ligand soluble, Rat, recombinant (sRANKL)	27
RANK Receptor soluble, Human, recombinant (sRANKR)	30
RANTES, Human, recombinant (CCL5)	71
RANTES, Mouse, recombinant (CCL5)	71
RANTES, Rat, recombinant (CCL5)	71
Relaxin-3, Human, recombinant	76
RELM - α , Mouse, recombinant	76
RELM - β , Human, recombinant	76
Resistin, Human, recombinant	77
Resistin, Mouse, recombinant	77
R-Spondin-1, Human, recombinant	64
S	
Semaphorin 3A, Human, recombinant	81
Stem Cell Factor, Human, recombinant (SCF)	21
Stem Cell Factor, Human, recombinant, Animal-derived-free (SCF)	21
Stem Cell Factor, Mouse, recombinant (SCF)	21
Stem Cell Factor, Mouse, recombinant, Animal-derived-free (SCF)	22
Stem Cell Factor, Rat, recombinant (SCF)	22
Stem Cell Growth Factor- α , Human, recombinant (SCGF- α)	22
Stromal Cell-Derived Factor-1 α , Human, recombinant (SDF-1 α /CXCL12)	68
Stromal Cell-Derived Factor-1 α , Human, recombinant, Animal-derived-free (SDF-1 α /CXCL12)	68

品名	ページ
Stromal Cell-Derived Factor -1 α , Mouse, recombinant (SDF-1 α / CXCL12)	68
Stromal Cell-Derived Factor -1 β , Human, recombinant (SDF-1 β / CXCL12)	68
Stromal Cell-Derived Factor -1 β , Mouse, recombinant (SDF-1 β / CXCL12)	68
T	
TARC, Human, recombinant (CCL17)	73
TECK, Human, recombinant (CCL25)	74
Thrombopoietin, Human, recombinant (TPO)	22
Thrombopoietin, Human, recombinant, Animal-derived-free (TPO)	22
Thrombopoietin, Mouse, recombinant (TPO)	22
Thrombopoietin, Rat, recombinant, Animal-derived-free (TPO)	23
TNF Receptor I soluble, Human, recombinant (sTNFR1)	30
TNF Receptor II soluble, Human, recombinant (sTNFR2)	30
TRAIL/Apo2L, Human, recombinant	28
TRAIL Receptor-1 soluble, Human, recombinant (sTRAILR1)	30
TRAIL Receptor-2 soluble, Human, recombinant (sTRAILR2)	31
Transforming Growth Factor - α , Human, recombinant (TGF- α)	37
Transforming Growth Factor - β 1, Human, recombinant (TGF- β 1)	60
Transforming Growth Factor - β 2, Human, recombinant(expressed in Insect Cells) (TGF- β 2)	61
Transforming Growth Factor - β 3, Human, recombinant (TGF- β 3)	61
Transforming Growth Factor - β 3, Human, recombinant, Animal-derived-free (TGF- β 3)	61
Trefoil Factor 1, Human, recombinant (TFF-1)	81
TSG, Human, recombinant	62
TSLP, Human, recombinant	23
Tumor Necrosis Factor - α , Human, recombinant (TNF- α)	27
Tumor Necrosis Factor - α , Human, recombinant, Animal-derived-free (TNF- α)	27
Tumor Necrosis Factor - α , Mouse, recombinant (TNF- α)	28
Tumor Necrosis Factor - α , Rat, recombinant (TNF- α)	28
Tumor Necrosis Factor - β , Human, recombinant (TNF- β)	28
TWEAK, Human, recombinant	29
TWEAK Receptor soluble, Human, recombinant (sTWEAKR)	31
V	
Vascular Endothelial Cell Growth Factor-A ₁₂₁ , Human, recombinant (VEGF-A ₁₂₁)	45
Vascular Endothelial Growth Factor-A ₁₂₁ , Human, recombinant, Animal-derived-free (VEGF-A ₁₂₁)	45
Vascular Endothelial Cell Growth Factor, Human, recombinant (VEGF-A ₁₆₅)	44
Vascular Endothelial Growth Factor-A ₁₆₅ , Human, recombinant, Animal-derived-free (VEGF-A ₁₆₅)	44
Vascular Endothelial Cell Growth Factor, Mouse, recombinant (VEGF-A ₁₆₅)	44
Vascular Endothelial Growth Factor-A ₁₆₅ , Rat, recombinant (VEGF-A ₁₆₅)	44
Vascular Endothelial Growth Factor-B ₁₆₇ , Human, recombinant (VEGF-B ₁₆₇)	45
Vascular Endothelial Growth Factor-C, Human, recombinant (VEGF-C)	45
Vascular Endothelial Growth Factor-D, Human, recombinant (VEGF-D)	46
VCAM-1, Human, recombinant	82
Visfatin, Human, recombinant	77
Vitronectin, Human, recombinant (VTN)	82
W	
Wnt-1, Human, recombinant	64
Wnt-3a, Mouse, recombinant, Solution	65
Wnt-7a, Human, recombinant	65



マウスES細胞の未分化能維持に

StemSure® LIF, マウス, 組換え体, 溶液

製品概要

●活性：10⁶units/mL

単位の定義：マウス ES 細胞株 (D3 株) を用いた細胞増殖促進アッセイにおいて、最大増殖度の50%の増殖度を与える量の1/20を1unitとする。

●推奨使用濃度：1,000units/mL

マウス ES 細胞株 (D3 株) の培養では、終濃度 1,000units/mL での使用を推奨

●形状：D-PBS, 1%BSA

●0.2µm フィルター滅菌済み

●保存条件：-20℃保存

※融解後は2～10℃保存の上、なるべく早くご使用下さい。




試験項目

- ◆ 実用試験 (マウス ES 細胞)
- ◆ 無菌試験
- ◆ エンドトキシン試験 など
- ◆ ALP 染色 (マウス ES 細胞)
- ◆ マイコプラズマ試験

データ

細胞形態・未分化マーカーの発現



1 × 10⁵ cells/10cm Dish で 5 日または 6 日後に継代を 14 回繰返し、位相差顕微鏡 (Phase) で撮影、ALP 染色、免疫染色 (各種未分化マーカー)、DAPI 染色を行った。

<培地組成>

D-MEM + 15%KSR + 4mmol/L L-Glutamine + 2×Non-essential Amino Acid + 1×Penicillin-Streptomycin + 0.1mmol/L StemSure® 2-Mercaptoethanol + 1,000units/mL StemSure® LIF

コードNo.	品名	規格	容量	希望納入価格(円)
199-16051	StemSure® LIF, Mouse, recombinant, Solution	細胞培養用	10 ⁶ units	30,000
195-16053			10 ⁶ units × 10	150,000



<http://www.wako-chem.co.jp/siyaku>

和光純薬工業株式会社

本 社：〒540-8605 大阪市中央区道修町三丁目1番2号 TEL: 06-6203-1788 (学術課)
支 店：〒103-0023 東京都中央区日本橋本町四丁目5番13号 TEL: 03-3270-8243 (学術課)
営業所：九州 TEL: 092-622-1005 中国 TEL: 082-285-6381 東海 TEL: 052-772-0788
筑波 TEL: 029-858-2278 東北 TEL: 022-222-3072 北海道 TEL: 011-271-0285

URL: <http://www.wako-chem.co.jp> E-mail: labchem-tec@wako-chem.co.jp
フリーダイヤル: 0120-052-099 フリーファックス: 0120-052-806

Wako Overseas Offices:

- Wako Chemicals USA, Inc.
Head Office (Richmond, VA) TEL: +1-804-714-1920 <http://www.wakousa.com>
Los Angeles Sales Offices (CA) TEL: +1-949-679-1700
Boston Sales Offices (MA) TEL: +1-617-354-6772
- Wako Chemicals GmbH (European Office) TEL: +49-2131-311-0 <http://www.wako-chemicals.de>